



## Sampling and preconcentration techniques for environmental radioactivity analysis

**Roos, Per**

*Publication date:*  
2009

[Link back to DTU Orbit](#)

*Citation (APA):*  
Roos, P. (Author). (2009). Sampling and preconcentration techniques for environmental radioactivity analysis. Sound/Visual production (digital)

---

### General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

# Sampling & preconcentration techniques for environmental radioactivity analysis



# Sampling techniques

- Air
- Soil
- Water
- Sediments





FIGURE 12.1—Degree of error in laboratory sample preparation relative to other activities


# When, where and how to sample ?

## An example from the dynamic Øresund area




Roos & Valeur, CSR,  
26(4), 474-487, 2006

**Continuous salinity (ppt)  
measurements at 9m depth stn 9  
(Nordre Röse)**



# Sampling strategy for Th-234 in water



## Rain-collector with ion-exchange column




Figure 2.8. Ion-exchange fallout collector.



Figure 2.7. Ion-exchange fallout collection and fallout pot collector.




## Wet & dry precipitation



Figure 2.9. Aerochem metrics wet/dry collector and EML wet/dry collector.



# Air – incomplete collection of all aerosols



Collection efficiency depends on

- Flow rate
- Total amount of collected material
- Filter type

Collecting resuspended soil & aerosols on sticky surfaces









# Soil sampling





# Radioactive particles

- Digital autoradiography, track-etch detectors and using dry or wet sample splitting using  $^{241}\text{Am}$  (59keV gamma).



# Techniques for sampling radioisotopes in water

- On site "bucket chemistry"
- Collecting water and bringing to home laboratory (only total conc !)
- In-situ sampling using sorbents




## Adsorbents for direct sampling in water

- Cs – Transition metal ferrocyanides.
- Th & Ra –  $\text{MnO}_2$
- Tc – Activated carbon





## Radiocesium

- Chloroplatinate ( $\text{H}_2\text{PtCl}_6$ )
- Ammoniummolybdophosphate
- Transition metal ferrocyanides



zink and copper ferrocyanide have been reported to have the highest caesium adsorption capacity reaching values of around 1.5 mol/mol for CuFC.

In normal environmental waters the concentrations of stable caesium is in the order of sub-ppb








Fig. 1. Large volume in-situ sampling system used for radiocesium. The equipment consists of one prefilter (1  $\mu\text{m}$ ), two CuFC impregnated, cotton-wound cartridge filters and a flow meter.



$$\text{Efficiency} = 1 - \frac{\text{CuFC(2)}}{\text{CuFC(1)}}$$




Fig. 2. The AMP /CuFC ratio as a function of filtered water volume for samples collected during the Swedish Antarctic Research




Fig. 3. The AMP/CuFC ratio as a function of ratio of radiocesium on the second and first CuFC filter,  $A_2/A_1$ . Samples from the Swedish Antarctic Research Expedition (SWEDARP), 1989.

# **MnO<sub>2</sub> impregnated fibre (sorption of Ra & Th)**



acrylic



polypropylene

## In-situ pumps for MnO<sub>2</sub> (or any other adsorbent)





# MnO<sub>2</sub> impregnated acrylic fibre



Gamma count

< 1 %

Melt/  
crush



5.3 %

ash

16.1 %

Gamma branching ratio:

63 keV  
3.8 %


93 keV  
5.4 %






automation

2004




AWI and ISITEC





# **Uptake of Th & Ra on $\text{MnO}_2$ limited by DOC in water**

**Volume of filtered water about 500 litres.**



## Tc-99 adsorption on activated carbon cartridges






Fig. 3. Adsorption efficiency for Tc from seawater on a cartridge active carbon filter.




# Sediment sampling




\*NIEMISTÖ, L. 1974. A gravity corer for studies of soft sediments. Merentutkimuslait. Julk./Havsforskningsinst. Skr. 238:33-38.



Core shortening (%) as a function of corer diameter (cm).  
From: Blomqvist, S. 1985. Reliability of core sampling of soft bottom sediment - an in situ study. Sedimentology (1985) 32, 605 - 612.




"Gemini" ready to descend. The two rudders and the valves are fixed with nylon strings and keep this position until bottom touch. A slight slack of the wire opens the pins and the rudders and valves are free to close.






The Box Corer is lowered towards the sea bed





The Box Corer is on its way up to the ship with the secured sample inside the tube



At the sea bed the corer tube penetrates into the sediment



The main cable is redrawn, and the closing grab goes into position beneath the corer tube



# Freeze-coring

166






Fig. 2.  $^{137}\text{Cs}$  in selected sediment sections and varve-counts for the two cores investigated. Fall-out from the Chernobyl accident resulted in a 7–10 fold increase of  $^{137}\text{Cs}$  in the Fe-rich layer deposited during the summer and autumn of 1986. A tentative correlation of individual varves between the two stations is indicated.



# Benthic landers



# Sediment traps and water samplers







## Samplers for speciation studies (chemical form, size etc)

### **Sediment & soil interstitial water**




- Centrifugation and filtering of soil or sediment interstitial water. Nitrogen flushing.
- Porous cups (with vacuum).
- Rhizon-samplers (artificial porous 'roots').
- Diffusive Gradient in Thin film samplers (DGT-samplers).

### **Water**


- Filtration & ultra-filtration.
- Selected co-precipitation techniques.
- Field Flow Fractionation (FFF)



# Rhizon-samplers



## Filtration 0.45 $\mu\text{m}$ with 293mm diameter membrane





# Cross-flow large volume ultrafiltration










# Ultrafiltration with permeate and retentate fraction







# In-situ filtration & ultrafiltration







**TABLE 10.1—Summary of sample preservation techniques.**

Preservation Technique	Advantages	Disadvantages
Addition of $\text{HNO}_3$	Reduces pH and inhibits plating of metals on container walls.	Strong oxidizer that might react with organic compounds, such as liquid scintillation cocktails. $^{14}\text{C}$ might be lost as $^{14}\text{CO}_2$ .
Addition of $\text{HCl}$	Reduces pH and inhibits plating of metals on container walls. Chloride forms strong anionic complexes with Iron and Uranium.	Causes quench in liquid scintillation cocktails. $^{14}\text{C}$ might be lost as $^{14}\text{CO}_2$ . Might cause corrosion of stainless steel planchets on gross analyses.
Addition of Sulfite	Forms a reducing environment to prevent the volatilization of iodine.	May produce undesirable oxidation states of iron or uranium.
Addition of Formaldehyde	Preserves organic samples. Prevents further biological activity.	May create disposal problems.
Cooling (Ice at approximately $0^\circ\text{C}$ )	Preserves organic samples (i.e., water, foods). Reduces dehydration and retains moisture. Reduces biological activity.	Ice melts, requiring replacement over time.
Freezing (Dry Ice at approximately $-78^\circ\text{C}$ )	Preserves organic samples (i.e., water, plant, animal). Suspends biological activity.	Dry ice sublimates and requires replacement. May crack sample container if frozen too quickly.
Addition of Paper Pulp	Provides large surface area for adsorption of metals, thus minimizing adsorption on container walls.	Requires pH to be one or less. Requires filtration and wet ashing of paper pulp and combining liquids to make a new solution.

# DGT-samplers



$$Q = \frac{gD \cdot A \cdot t}{\Delta g} \cdot C_b$$





# Ship-based Field Flow Fractionation system



# Sample treatment

- Pre-treatment involving drying, sieving, ashing, filtering and precipitation
- Chemical preparation, purification and isolation of the elements in question



# Measurement methods & radiochemistry

	Chemical analysis	Liquid scintillation counting	Gamma spectrometry	Total beta counting	Alpha spectrometry	ICPMS	Neutron activation analysis
$^3\text{H}$	X	X					
$^{90}\text{Sr}$	X	X		X			
$^{99}\text{Tc}$	X			X		X	
$^{129}\text{I}$	X						X
$^{137}\text{Cs}$	X		X				
$^{210}\text{Pb}$			X				
$^{210}\text{Po}$	X				X		
$^{226}\text{Ra}$	X	X	X				
$^{230}\text{Th}$	X				X		
$^{232}\text{Th}$	X				X	X	
$^{234}\text{U}$	X				X	X	
$^{238}\text{U}$	X				X	X	
$^{237}\text{Np}$	X				X	X	
$^{239}\text{Pu}$	X					X	
$^{240}\text{Pu}$	X					X	
$^{239+240}\text{Pu}$	X	X			X		
$^{241}\text{Am}$	X				X		
Other gamma emitters			X				



# Decomposition methods



07-maj-08 17:17



- **Dry oxidation** -550-600 °C. Too high temperature-volatile elements, oxides difficult to dissolve.
- Use of ashing acids,  $\text{HNO}_3$ ,  $\text{H}_2\text{SO}_4$ , oxalic acid

**TABLE 12.3 — Examples of dry-ashing temperatures (platinum container)**

Element	Temperature/Matrix
Cobalt	450–600 °C for biological material; some losses reported due to reactions with crucible; increased volume of sample increases volume of ash and limits loss of sample.
Cesium	400–450 °C for food and biological material; CsCl and CsNO <sub>3</sub> begin to volatilize when held at temperatures above 500 °C for any length of time.
Iodine	450–500 °C with an alkaline ashing aid to prevent volatilization; losses reported for temperatures as low as 450 °C even with alkaline ashing aids added; total volatilization >600 °C.
Lead	450–500 °C acceptable for most samples; bone or coal (lead phosphate) may be ashed as high as 900 °C without significant losses; PbO <sub>2</sub> reacts with silica in porcelain glaze at low temperatures; PbCl <sub>2</sub> is relatively volatile and nitrate or sulfate ashing aids have been used to good effect.
Plutonium	450 °C with nitric acid ashing aid for biological material, 550 °C for dust on air filters, 700 °C for soil; high temperature leads to adsorption onto carbon particles and incomplete dissolution of ash.
Strontium	450–550 °C for plants, 600 °C for meat, 700 °C for milk and bone.
Technetium	725–750 °C for plants treated with ammonia.
Thorium	750 °C for bone.
Uranium	600 °C for coal, 750 °C for biological material; uranium reacts with porcelain glaze resulting in sample losses.

Source: Bock (1979).

(Note that reducing conditions for platinum containers are given in Table 12.2)

**TABLE 12.1 — Examples of volatile radionuclides**

Isotope	Chemical Form	Boiling Point (°C) *
Tritium — $^3\text{H}$	$\text{H}_2\text{O}$	$100^\circ$
Carbon — $^{14}\text{C}$	$\text{CO}_2$ (produced from $\text{CO}_3^{-2}$ or oxidation of organic material)	$-78.5^\circ$
	Magnesium, calcium, and sodium carbonates	Natural ores of these metals decompose between $825^\circ$ and $1,330^\circ$ to yield the respective metal oxides
Iodine — $^{131}\text{I}$ , $^{129}\text{I}$	$\text{I}_2$	$185.2^\circ$ (sublimes readily)
Cesium — $^{134}\text{Cs}$ , $^{135}\text{Cs}$ , $^{136}\text{Cs}$ , $^{137}\text{Cs}$	$\text{Cs}^0$ (as metal)	$678.4^\circ$ (melts at $28^\circ$ )
	$\text{Cs}_2\text{O}$ (as metallic oxide) (nitrates decompose to oxides)	$\sim 400^\circ$
	$\text{CsCl}$ (as metallic chloride)	$1290^\circ$
Technetium — $^{99}\text{Tc}$	$\text{Tc}_2\text{O}_7$	$310.6^\circ$
	$\text{TcCl}_4$	Sublimes above $300^\circ$
	$\text{TcO}_2$	Sublimes above $900^\circ$
	[Most Tc compounds sublime above $300^\circ$ . Tc(VII) is an oxidant that reacts with organic solvents forming Tc(IV)]	
Polonium — $^{208}\text{Po}$ , $^{209}\text{Po}$ , $^{210}\text{Po}$	$\text{Po}^0$	$962^\circ$
	$\text{PoCl}_4$	$390^\circ$
	$\text{Po}(\text{NO}_3)_4$ [as a solid]	Decomposes to $\text{PoO}_2$ above $\sim 150^\circ$
	$\text{PoO}_2$	Decomposes to Po metal above $500^\circ$
Lead — $^{210}\text{Pb}$ , $^{212}\text{Pb}$ , $^{205}\text{Pb}$	$\text{Pb}^0$	$1744^\circ$
	$\text{PbCl}_2$	$950^\circ$
	$\text{Pb}(\text{NO}_3)_2$	Decomposes to oxide above $470^\circ$
	$\text{PbO}$	$888^\circ$

\* The closer the sample preparation temperature is to the boiling point of the compound, the more significant will be the loss of the material. However, if the objective is to distill the analyte compound from other nonvolatile materials, then boiling temperature is needed. Sample preparation near the decomposition temperature should be avoided for those compounds that have a decomposition temperature listed in the table.

Sources: Greenwood and Earnshaw (1984); Windholz (1976); Schwochau (2000); Sneed and Brasted (1958).



## Wet oxidation systems

- $\text{HNO}_3\text{-HClO}_4$ —danger of explosions
- $\text{HNO}_3\text{-H}_2\text{SO}_4$
- $\text{HNO}_3\text{-HCl-H}_2\text{O}_2$
- The last two ones are safer but do not decompose all organic compounds
- **Oxidative fusion**
- For minerals, sediment, fiber glass filters
- $\text{Na}_2\text{O}_2$ ,  $\text{NaOH}$ ,  $\text{KHSO}_4$ ,  $\text{K}_2\text{S}_2\text{O}_7$
- Often very violent reactions

## **Wet oxidation**

- **-Advantages:**
- lower temperature
- improves isotopic exchange
- low pH,  $H^+$  ions fill vacancies in wall containers
- **Disadvantages**
- formation of compounds difficult to dissolve
- background from radionuclides in acids and containers

**TABLE 13.2 — Examples of acids used for wet ashing**

Acid	Typical Uses
Hydrofluoric Acid, HF	Removal of silicon and destruction of silicates; dissolves oxides of Nb, Ta, Ti, and Zr, and Nb, and Ta ores.
Hydrochloric Acid, HCl	Dissolves many carbonates, oxides, hydroxides, phosphates, borates, and sulfides; dissolves cement.
Hydrobromic Acid, HBr	Distillation of bromides (e.g., As, Sb, Sn, Se).
Hydroiodic Acid, HI	Effective reducing agent; dissolves $\text{Sn}^{+4}$ oxide and $\text{Hg}^{+2}$ sulfide.
Sulfuric Acid, $\text{H}_2\text{SO}_4$	Dissolves oxides, hydroxides, carbonates, and various sulfide ores; hot concentrated acid will oxidize most organic compounds.
Phosphoric Acid, $\text{H}_3\text{PO}_4$	Dissolves $\text{Al}_2\text{O}_3$ , chrome ores, iron oxide ores, and slag.
Nitric Acid, $\text{HNO}_3$	Oxidizes many metals and alloys to soluble nitrates; organic material oxidized slowly.
Perchloric Acid, $\text{HClO}_4$	Extremely strong oxidizer; reacts violently or explosively to oxidize organic compounds; attacks nearly all metals.



**TABLE 13.1 — Common fusion fluxes**

Flux (mp, °C)	Fusion Temperature, °C	Type of Crucible	Types of Sample Decomposed
$\text{Na}_2\text{S}_2\text{O}_7$ (403°) or $\text{K}_2\text{S}_2\text{O}_7$ (419°)	Up to red heat	Pt, quartz, porcelain	For insoluble oxides and oxide-containing samples, particularly those of Al, Be, Ta, Ti, Zr, Pu, and the rare earths.
NaOH (321°) or KOH (404°)	450-600°	Ni, Ag, glassy carbon	For silicates, oxides, phosphates, and fluorides.
$\text{Na}_2\text{CO}_3$ (853) or $\text{K}_2\text{CO}_3$ (903)	900-1,000°	Ni Pt for short periods (use lid)	For silicates and silica-containing samples (clays, minerals, rocks, glasses), refractory oxides, quartz, and insoluble phosphates and sulfates.
$\text{Na}_2\text{O}_2$	600°	Ni; Ag, Au, Zr; Pt (<500 °C)	For sulfides; acid-insoluble alloys of Fe, Ni, Cr, Mo, W, and Li; Pt alloys; Cr, Sn, and Zn minerals.
$\text{H}_3\text{BO}_3$	250°	Pt	For analysis of sand, aluminum silicates, titanite, natural aluminum oxide (corundum), and enamels.
$\text{Na}_2\text{B}_4\text{O}_7$ (878°)	1,000-1,200°	Pt	For $\text{Al}_2\text{O}_3$ , $\text{ZrO}_2$ and zirconium ores, minerals of the rare earths, Ti, Nb, and Ta, aluminum-containing materials; iron ores and slags.
$\text{Li}_2\text{B}_4\text{O}_7$ (920°) or $\text{LiBO}_2$ (845°)	1,000-1,100°	Pt, graphite	For almost anything except metals and sulfides. The tetraborate salt is especially good for basic oxides and some resistant silicates. The metaborate is better suited for dissolving acidic oxides such as silica and $\text{TiO}_2$ and nearly all minerals.
$\text{NH}_4\text{HF}_2$ (125°) NaF (992°) KF (857°) or KHF <sub>2</sub> (239°)	900°	Pt	For the removal of silicon, the destruction of silicates and rare earth minerals, and the analysis of oxides of Nb, Ta, Ti, and Zr.

Source: Dean (1995) and Bock (1979).

# Preconcentration techniques

- Volatilisation & distillation (I, Ru, Po, Ni, Tc etc) – limited nr of elements
- Precipitation and co-precipitation techniques.

**TABLE 14.13 — Common coprecipitating agents for radionuclides<sup>(1)</sup>**

Radionuclide	Oxidation State	Coprecipitate	Carrier <sup>(2)</sup>	Notes
Am	+3	hydroxide iodate fluoride, oxalate, phosphate, hydroxide oxalate acetate fluoride, sulfate acetate	Am <sup>+3</sup> , Fe <sup>+3</sup> Ce <sup>+4</sup> , Th <sup>+4</sup> , Zr <sup>+4</sup> La <sup>+3</sup> , Ce <sup>+3</sup> , Nd <sup>+3</sup> , Bi <sup>+3</sup>  Ca <sup>+2</sup> Am <sup>+4</sup> La <sup>+3</sup> UO <sub>2</sub> <sup>+2</sup>	
Cs	+1	phosphomolybdate, chloroplatinate, bismuth nitrate, silicomolybdate	Cs <sup>+1</sup>	
Co	+2	hydroxide potassium cobalt nitrate 1-nitroso-2-naphthol sulfide	Co <sup>+2</sup> Co <sup>+2</sup> Co <sup>+2</sup> Co <sup>+2</sup>	
Fe	+3	hydroxide ammonium pyrouanate	Fe <sup>+3</sup> Fe <sup>+3</sup>	
I	-1	iodide	Pb <sup>+2</sup> , Ag <sup>+1</sup> , Pd <sup>+2</sup> , Cu <sup>+2</sup>	
Ni	+2	dimethylglyoxime hydroxide	Ni <sup>+2</sup>	
Nb	(V)	hydroxide, phosphate	Nb(V)	
Np	+4	phosphate	Ca <sup>+2</sup>	
Po	+4	tellurium tellurate selenium dioxide hydroxide sulfide	Te Pb <sup>+2</sup> Se or Se <sup>-2</sup> Mn <sup>+4</sup> Fe <sup>+3</sup> , Al <sup>+3</sup> , La <sup>+3</sup> Cu <sup>+2</sup> , Bi <sup>+2</sup> , Pb <sup>+2</sup>	Tellurate reduced with SnCl <sub>2</sub>
Pu	+3 +4 (VI)	fluoride sulfate fluoride oxalate, iodate phosphate sodium uranylacetate	La <sup>+3</sup> , Nd <sup>+3</sup> , Ce <sup>+3</sup> , Ca <sup>+2</sup> La <sup>+3</sup> (K <sup>+1</sup> ) La <sup>+3</sup> , Nd <sup>+3</sup> , Ce <sup>+3</sup> Th <sup>+4</sup> Zr <sup>+2</sup> , Bi <sup>+3</sup> UO <sub>2</sub> <sup>+2</sup>	



Radionuclide	Oxidation State	Coprecipitate	Carrier <sup>(2)</sup>	Notes
Ra	+2	hydroxide sulfate, chromate, chloride, bromide oxalate, phosphate fluoride	Fe <sup>+3</sup> Ba <sup>+2</sup> Th <sup>+4</sup> , Ca <sup>+2</sup> , Ba <sup>+2</sup> La <sup>+3</sup>	
Sr	+2	carbonate nitrate chromate sulfate phosphate hydroxide	Sr <sup>+2</sup> , Ba <sup>+2</sup> , Ca <sup>+2</sup> Sr <sup>+2</sup> , Ba <sup>+2</sup> Ba <sup>+2</sup> Sr <sup>+2</sup> , Ca <sup>+2</sup> , Pb <sup>+2</sup> Sr <sup>+2</sup> Fe <sup>+3</sup>	Alkaline pH
Tc	+4 (VII)	hydroxide chlorate, iodate, perruthenate, tetrafluoroborate sulfide	Tc <sup>+4</sup> , Fe <sup>+3</sup> , Mn <sup>+2</sup> (Phenyl) <sub>4</sub> As <sup>+1</sup> Tc <sup>+7</sup> , Re <sup>+7</sup> , Cu <sup>+2</sup> , Cd <sup>+2</sup>	
Th	+4	hydroxide fluoride iodate phosphate, peroxide sulfate oxalate	Th <sup>+4</sup> , La <sup>+3</sup> , Fe <sup>+3</sup> , Zr <sup>+3</sup> , Ac <sup>+3</sup> , Zn <sup>+2</sup> Th <sup>+4</sup> , La <sup>+3</sup> , Nd <sup>+3</sup> , Ce <sup>+3</sup> Th <sup>+4</sup> , Zr <sup>+3</sup> Th <sup>+4</sup> , Bi <sup>+3</sup> Ba <sup>+2</sup> Ca <sup>+2</sup>	
U	+4	cupferron, pyrophosphate, phosphate, iodate, sulfate, oxalate	U <sup>+4</sup>	
		fluoride	La <sup>+3</sup> , Nd <sup>+3</sup>	
	(V)	phosphate	Zr <sup>+3</sup>	
		sulfate	Ca <sup>+2</sup>	
	(VI)	cupferron	U(VI)	Neutral solution
		pyrouranate	U(VI)	From aqueous NH <sub>3</sub> , many ions stay in solution as NH <sub>3</sub> complex
		phosphate	U(VI), Al <sup>+3</sup>	
		peroxide	U(VI)	Th <sup>+4</sup> , Zr <sup>+3</sup> also coprecipitate
		hydroxide	Fe <sup>+3</sup>	Without carbonate
		fluoride	Th <sup>+4</sup>	
Zr	+4	hydroxide	Fe <sup>+3</sup>	

(1) Compiled from: Anders, 1960; Booman and Rein, 1962; Cobble, 1964; EPA, 1973; 1980; 1984; DOE, 1990, 1995, 1997; Finston and Kinsley, 1961, Grimaldi, 1961; Grindler, 1962; Hyde, 1960; Kallmann, 1961; Kallmann, 1964; Kirby and Salutsky, 1964; Metz and Waterbury, 1962; Sedlet, 1964; Sundermann and Townley, 1960; and Turekian and Bolter, 1966.

(2) If the radionuclide itself is listed as the carrier, a different isotope would be used to assess recovery.

**TABLE 14.14 — Coprecipitation behavior of plutonium and neptunium**

Carrier Compound	Pu <sup>+3</sup>	Pu <sup>+4</sup>	Pu(VI)	Np <sup>+4</sup>	Np(V)	Np(VI)
Hydroxides	C	C	C	C	C	C
Calcium fluoride	C	C		C		
Lanthanum fluoride	C	C	NC	C	C	NC
Barium sulfate	C	C	NC	C	NC	NC
Phosphates:						
Calcium phosphate	C	C		C		
Bismuth phosphate	C	C		C	NC	NC
Zirconium phosphate	NC	C	NC	C	NC	NC
Thorium pyrophosphate	NC	C	NC			
Thorium hypophosphate		C	NC			
U <sup>+4</sup> hypophosphate		C	NC			
Oxalates:						
Lanthanum oxalate	C	C	NC	NC		
Bismuth oxalate	C	C	NC			
Thorium oxalate	C	C	NC	C		
U <sup>+4</sup> oxalate	C	C	NC			
Iodates:						
Zirconium iodate		C	NC	C		
Ceric iodate		C	NC	C		
Thorium iodate		C	NC	C		NC
Sodium uranyl acetate	NC	NC	C	NC	Poor	C
Zirconium phenylarsenate	NC	C	NC	C	Poor	NC
Thorium peroxide		C		C		
Bismuth arsenate		C	NC	C		


“C” indicates nearly quantitative coprecipitation under proper conditions; “NC” indicates that coprecipitation can be made less than 1–2 percent under proper conditions. [Data compiled from Seaborg and Katz, Korkisch (1969), and the NAS-NS 3050, 3058 and 3060 monographs.]

# MARLAP-manual from US EPA

NUREG-1576  
EPA 402-B-04-001B  
NTIS PB2004-105421



## Multi-Agency Radiological Laboratory Analytical Protocols Manual Volume II: Chapters 10 – 17 and Appendix F



July 2004

<http://www.epa.gov/rpdweb00/marlap/manual.html>