



## **Growth Dynamics of Layered Carbon Studied by In Situ Environmental Transmission Electron Microscopy**

**Kling, Jens**

*Publication date:*  
2015

*Document Version*  
Peer reviewed version

[Link back to DTU Orbit](#)

*Citation (APA):*

Kling, J. (2015). *Growth Dynamics of Layered Carbon Studied by In Situ Environmental Transmission Electron Microscopy*. Abstract from 2015 MRS Fall Meeting and Exhibit, Boston, Massachusetts, United States.

---

### **General rights**

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

## Growth dynamics of layered carbon studied by *in situ* environmental transmission electron microscopy

Nanostructured carbon materials are predicted to play a major role in future electronic applications. Cheaper and smaller components with improved or new functionality and lower power consumption are necessary, where conventional materials reach their limitations. Devices based on layered carbon materials, such as graphene or multilayer graphene, open for a giant technological leap in electro and optical performance. A cheap way to synthesize such materials on a large scale is chemical vapor deposition (CVD) growth on catalyst films like copper or nickel. However, the growth processes and the control of these is still not completely understood.

We present *in situ* transmission electron microscopy (TEM) experiments in a FEI Titan 80-300 Environmental TEM (ETEM) for studying the growth of layered carbon materials on these catalysts. The ETEM allows imaging with controlled gas environments around the sample up to a few mbar. In combination with a MEMS-based heating holder, growth of layered carbon materials is systematically studied at the atomic level.

NiO particles in the size range up to a few hundred nm are used as model material for the large-scale growth, as they provide a small surface curvature and mimic the native oxidized Ni surface found in conventional Ni thin films. The particles are reduced in the microscope under H<sub>2</sub> at 500-600°C in order to form a catalytically active Ni surface. Introducing carbon precursor gas at elevated temperature leads to growth of layered carbon. By following the appearance of carbon layers, the layer growth rate as well as the growth rate of individual graphene layers on the Ni surface can be determined directly from the ETEM observations.

These results add to a better understand the growth mechanism and help to control and optimize the formation process.