



## **Speciation analysis of trace elements in food and feed - status and future developments**

**Sloth, Jens Jørgen**

*Publication date:*  
2011

*Document Version*  
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

*Citation (APA):*  
Sloth, J. J. (2011). *Speciation analysis of trace elements in food and feed - status and future developments*. Poster session presented at AOAC 125th Annual Meeting & Exposition, New Orleans, USA.

---

### **General rights**

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal


If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

# Speciation analysis of trace elements in food and feed - status and future developments

---

*Jens J. Sloth*

National Food Institute (DTU Food)  
Technical University of Denmark



# AGENDA

---

- Speciation in relation to food and feed safety
- Selected examples
  - 1. Arsenic speciation analysis (importance of inorganic arsenic)
  - 2. Organotin speciation analysis (food contact materials)
  - 3. Selenium speciation analysis (food fraud!)
  - .....and a little bit of nano!



# CURRENT SITUATION IN EU LEGISLATION:


## Foodstuffs

MLs for Pb, Cd, Hg and Sn  
EU directive 2006/1881/EC

## Animal feedingstuffs


MLs for As, Pb, Cd and Hg  
EU directive 2002/32/EC

Only maximum levels for  
TOTAL CONCENTRATION  
of the metals



## Example – arsenic speciation

### Important for correct risk assessment



~ 40 – 250  $\mu\text{g}$  inorg As

1 kg rice => 50-300  $\mu\text{g}$  As



<10  $\mu\text{g}$  inorg As

1 kg fish => 3000-10000  $\mu\text{g}$  As

There is most focus on rice from a food safety point of view – why???



The chemical form of arsenic is important  
=> Arsenic speciation

Focus  
on  
Food  
Safety

# EFSA (2009) and JECFA (2010) opinions on arsenic in food

- Old PTWI value (WHO, 1988) was withdrawn ( $15 \mu\text{g/kg bw/week}$ )




- **NEW!**  $\text{BMDL}_{1.0} = \underline{0.3 - 8 \mu\text{g/kg bw per day}}$  for inorganic arsenic
- $\Rightarrow$  EU dietary exposures within this range for average and high level consumers
- $\Rightarrow$  Risk to some consumers cannot be excluded



- **NEW!**  $\text{BMDL}_{0.5} = \underline{3 \mu\text{g/kg bw per day}}$  for inorganic arsenic
- $\Rightarrow$  0.5% increased incidence of lung cancer for 12 y exposure

- “...there is a need to produce speciation data for different food commodities to support dietary exposure assessment...”
- “...more accurate information on the inorganic arsenic content of foods is needed to improve assessments of dietary exposures to inorganic arsenic”
- “...need for validated methods for selective determination of inorganic arsenic in food matrices”

# Inorganic arsenic in wild caught fish => no concern



## Norwegian survey

900 individual fish samples

- Atlantic halibut
- Cod
- Greenland halibut
- Mackerel
- Herring
- Tusk

## Results

Total arsenic.....0.3-110 mg/kg


Inorganic arsenic.... < 0.01 mg/kg

(only 37 samples > LOQ)



N I F E S  
NATIONAL INSTITUTE  
OF NUTRITION AND  
SEAFOOD RESEARCH

...but in bivalves high contents in some samples...




Data from 175 blue mussel (*Mytilus edulis*) samples collected along the Norwegian Coastline.





# Arsenic in rice products DK - survey



## 105 samples

- white rice
- brown rice
- red rice
- black rice
- rice crackers


## 33 samples > 0,15 mg/kg

- 2 parboiled (20%)
- 4 brown (50%)
- 4 red (50%)
- 5 black (71%)
- 1 Basmati (10%)
- 1 Pudding rice (9%)
- 1 wild rice (20%)
- 15 rice crackers (100%)

# Arsenic in rice products DK - survey



iAs vs total As



## 105 samples

- white rice
- brown rice
- red rice
- black rice
- rice crackers

## 22 samples > 0,2 mg/kg

- 1 parboiled (10%)
- 1 brown (12%)
- 3 red (37%)
- 2 black (28%)
- 0 Basmati (0%)
- 0 Pudding rice (0%)
- 0 wild rice (0%)
- 15 rice crackers

Rice cracker mean: 0.31 mg/kg – intake 50 g/dag => 15 µg iAs (~1 µg/kg bw/dag)  
> EFSA BMDL<sub>01</sub> 0.3-8 µg/kg bw/dag



# Inorganic arsenic in chinese food supplements


Name of Food supplement	Total Arsenic ( $\mu\text{g/g}$ )	Inorganic arsenic ( $\mu\text{g/g}$ )
Xiao Yao Wan	0.82	0.85
Bu Zhong Yi Qi Wan	0.62	0.50
Da Bu Yin Wan	0.59	0.55
Six Flavor teapills	0.72	N.D.
Golden Book Teapills	0.58	0.57
Xiang Sha Liu Jun Zi Wan	0.94	0.80
<b>Gan Mao Ling</b>	<b>1.24</b>	<b>1.01</b>
Chuan Xin Lian	5.00	3.17
Bi Yan Pian	0.70	0.58
Arouse power	1.12	1.02
Bio Chlorella	0.62	0.21
Unik Spirulina Kapsler	2.59	0.13
Chlorella	0.58	0.03
Ez-Biloba	0.63	0.67
Qvinde Dong Quai	0.68	0.48



Gan Mao Ling  
(against flu and common cold)

Rec dose: 18 pills per day  
 $\Rightarrow$  iAs  $\sim 13 \mu\text{g/day}$   
 $\Rightarrow 0.22 \mu\text{g/kg bw/day}$  (@60 kg)  
 Close to EFSA  $\text{BMDL}_{01}$  !!

# Inorganic As in SRM NIST1586a – literature data



Extraction solution	Detection	iAs (µg/kg)	Reference
1 2M TFA	LC-ICPMS	87 +/- 9	Ackermann (2005)
2 Enzymatic, pepsin and pancreatin	LC-ICPMS	101 +/- 7	Ackermann (2005)
3 2M TFA	LC-ICPMS	92 +/- 2	Heitkemper (2009)
4 2M TFA	LC-ICPMS	80 +/- 16	Williamson (2008)
5 2M TFA	LC-ICPMS	100 +/- 10	Williamson (2008)
6 MeOH:H2O	LC-ICPMS	109 +/- 3	D'Amato (2008)
7 Enzymatic, alfa-amylase	LC-ICPMS	106 +/- 7	Kohlmeier (2009)
8 Enzymatic, protease and alfa-amylase	LC-ICPMS	88 +/- 6	Sanz (2009)
9 1M H3PO4	HG-AFS	102 +/- 2	Matosreyes (2008)
10 1% HNO3	LC-ICPMS	99 +/- 4	Raab (2009)
11 1% HNO3	LC-ICPMS	99 +/- -	Sun (2008)
12 1% HNO3	LC-ICPMS	110 +/- 10	Sun (2009)
13 0,5 M TFA	LC-ICPMS	77 +/- -	Heitkemper (2009)
14 enzymatic, proteas and alfa-amylase	LC-ICPMS	96 +/- 9	Mar (2009)
15 2M TFA	LC-ICPMS	100 +/- 12	Meharg (2008)
16 water	LC-ICPMS	96 +/- 3	Narukawa (2008)
17 0,07M HCL and 10%H2O2	LC-ICPMS	103 +/- 15	DTU Food (2009)

Good agreement  
between  
labs and methods

# Animal feedingstuffs

## Commission directive 2009/114/EC (amendment)

Undesirable substances	Products intended for animal feed	Maximum content in mg/kg (ppm) relative to a feedingstuff with a moisture content of 12 %
(1)	(2)	(3)
1. Arsenic (*) (**)	Feed materials with the exception of:	2
	— meal made from grass, from dried lucerne and from dried clover, and dried sugar beet pulp and dried molasses sugar beet pulp,	4
	— palm kernel expeller,	4 (***)
	— phosphates and calcareous marine algae,	10
	— calcium carbonate,	15
	— magnesium oxide,	20
	— feedingstuffs obtained from the processing of fish or other marine animals, including fish,	25 (***)
	— seaweed meal and feed materials derived from seaweed,	40 (***)
	Iron particles used as tracer.	50
	Additives belonging to the functional group of compounds of trace elements except:	30



Only max levels for total arsenic

### FOOTNOTE

(\*\*\*) Upon request of the competent authorities, the responsible operator must perform an analysis to demonstrate that the content of inorganic arsenic is lower than 2 ppm in feedingstuffs of particular importance for the seaweed species *Hizikia fusiforme*.

— zinc oxide, manganese oxide and copper oxide,	100
Complete feedingstuffs with the exception of:	2
— complete feedingstuffs for fish and complete feedingstuffs for fur animals,	10 (***)
Complementary feedingstuffs with the exception of:	4
— mineral feedingstuffs,	12


Speciation analysis is required !!

# SPE-HG-AAS — a speciation alternative...

μ-wave  
extraction

Separation  
by SPE

Detection by  
HG-AAS



Sequential elution for selective off-line separation of inorg As from organo As species by SPE

Abundance

75000  
70000  
65000  
60000  
55000  
50000  
45000  
40000  
35000  
30000

Sample  
load

OrganoAs

Wash

HPLC-ICPMS of SPE fractions

Sample eluate  
inorganic As



# SPE-HG-AAS

μ-wave extraction

Separation by  
SPE

Detection by HG-  
AAS

## In-house validation data



Parameter	Result
Analysis time	2 x 7 h for 24 samples
LoD (mg/kg)	0.08
LoQ (mg/kg)	0.16
Repeatability (%RSD)	3 - 7
Accuracy (%)	90 - 104

## Collaborative trial in marine feed material



- 10 labs provided results
- Method working range tested: 0.1 – 2.6 mg Kg<sup>-1</sup>
- HorRat values <2 in the working range tested
- The method is fit for purpose
- Final publication as EN expected July 2012




# Routes of exposure to organotin compounds

Used in

- Agriculture
- Antifoulings
- Industry

- PVC-Materials

**TDI: 0.25  $\mu\text{g}/\text{kg bw}/\text{day}$**   
 **$\sum$  TBT, DBT, TPhT and DOT**







# Legislation on OTCs in Food Contact Materials


Compounds	Maximum level ( $\mu\text{g Sn/kg foodstuff}$ )
$\sum$ DBT, TBT, TPhT and DOT	40 (6)
$\sum$ MMT, DMT	180
MOT	1200
MDDT	12000 (50)
DDDT	24000 (50)

Ref: EFSA (2005); proposed EFSA values in parenthesis

## Assumptions:

- 1 kg food per 6 dm<sup>2</sup>
- 100 mL in contact with 0.6 dm<sup>2</sup>

- Max levels on organotin migrating from the packaging material
- Testing by the use of food simulators (water, acid, oil, alcohol etc)
- BUT no maximum levels on organotins in the foodstuff itself!!




# Organotin migration from Food Contact Materials II

## Small scale survey on 33 FCMs

Baking paper, PVC cling films, silicone baking forms, lids with PVC gaskets


PUR-agglomerated cork wine stoppers



- PVC lid
- 3% acetic acid

## Overlaid standard and sample

➤ DBT concentration: 9.9 µg/kg



Output of DK survey:

- 33 samples
- 11 contained OT (mainly DBT)
- 3 exceeded EFSA guideline limit


# Selenium in commercial food supplements

*Organic bound Se?*



*Selenite?*

*Selenate?*



*Se yeast?*


*Selenomethionine?*

*Amino chelated Se?*

**Is the selenium source declared correctly??**

# Se speciation by HPLC-ICPMS

Supplement declared as organic bound Selenium (125 µg/tablet)



# Speciation summary

---



- ✓ speciation methods are more and more commonly used
- ✓ instrumentation is widely available
- ✓ legislation on species has started
- ✓ ...and more is expected in the future!



- ✓ standardised methods are not ready!
- ✓ ...but the need is known by authorities
- ✓ legislation is still behind!
- ✓ Lack of CRMs (e.g. for iAs)

# Perspective – trace elements in nanoform

## NanoLyse Project "Nanoparticles in Food: Analytical methods for detection and characterisation"

Validated methods for the determination of inorganic ENP in food extracts, based on size separation, size determination and specific detection




*silver nanoparticles in lean meat*



*silica nanoparticles in tomato soup*



# DTU Food Analytical platform: AF4-MALS/DLS-ICP-MS



asymmetric flow  
field flow  
fractionation  
(AF4)

optical detection  
multi angle (MALS)  
and dynamic light  
scattering (DLS), UV-  
vis absorption and  
fluorescence


inductively coupled  
plasma mass  
spectrometry  
(ICP-MS)

particle separation  
according to their  
size (1nm – few  $\mu$ m)

particle detection  
size determination

elemental detection  
chemical identity  
quantification

Fractogram =>



## Acknowledgements and funding sources

### Coworkers:

Rie R. Rasmussen, Rikke V. Hedegaard, Bjørn Schmidt,  
Xenia T. Trier, Katrin Löschner and Erik H. Larsen



Inge Rokkjær, Gudrun Hilbert and Dorte L. Cederberg



Kåre Julshamn and A.K. Lundebye



Fernando Correiro-Raposo, Håkan Emteborg and Beatriz de la Calle



### Funding sources:

European Community's Seventh Framework Programme



European Committee for standardisation (CEN)



Danish Food Administration (DFVA)



125<sup>th</sup> AOAC

ANNUAL MEETING & EXPOSITION



Thanks for your attention! 