



## **Automation and Methodology Development for Environmental and Biological Determination of Pu, Np, U and Tc**

**Qiao, Jixin**

*Publication date:*  
2013

[Link back to DTU Orbit](#)

*Citation (APA):*

Qiao, J. (Author). (2013). Automation and Methodology Development for Environmental and Biological Determination of Pu, Np, U and Tc. Sound/Visual production (digital)  
[http://www.nks.org/en/seminars/presentations/nks-b\\_radioanalysis\\_workshop\\_2-6\\_september\\_2013.htm](http://www.nks.org/en/seminars/presentations/nks-b_radioanalysis_workshop_2-6_september_2013.htm)

---

### **General rights**

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

# **Automation and Methodology Development for Environmental and Biological Determination of Pu, Np, U and Tc**

Jixin Qiao

Radioecology and Tracer Studies  
DTU Nutech

06-09-2013

# Properties of Pu, Np, U and Tc

Nuclide	Isotope	Main origination	Half-life	Main production	Principal decay mode
Pu	$^{238}\text{Pu}$	Anthropogenic	87.7 y	NA and $\beta$ decay of $^{235}\text{U}$ and $^{238}\text{U}$	$\alpha$
	$^{239}\text{Pu}$		$2.4 \times 10^4$ y	Bombardment of $^{238}\text{U}$	$\alpha$
	$^{240}\text{Pu}$		$6.6 \times 10^3$ y	$^{239}\text{Pu}$ (n, $\gamma$ ) $^{240}\text{Pu}$	$\alpha$
	$^{241}\text{Pu}$		14.4 y	$^{240}\text{Pu}$ (n, $\gamma$ ) $^{241}\text{Pu}$	$\beta^-$
Np	$^{237}\text{Np}$	Anthropogenic	$2.4 \times 10^6$ y	NA and $\beta$ decay of $^{235}\text{U}$ and $^{238}\text{U}$	$\alpha$
U	$^{234}\text{U}$	Natural	$2.4 \times 10^6$ y		$\alpha$
	$^{235}\text{U}$	Natural	$2.5 \times 10^5$ y		$\alpha$
	$^{236}\text{U}$	Anthropogenic	$2.3 \times 10^7$ y	$^{235}\text{U}$ neutron activation (NA)	$\alpha$
	$^{238}\text{U}$	Natural	$4.5 \times 10^9$ y		$\alpha$
Tc	$^{99}\text{Tc}$	Anthropogenic	$2.1 \times 10^5$ y	$^{235}\text{U}$ , $^{239}\text{Pu}$ fission product	$\beta^-$

# Sources of Pu, Np, U and Tc in the environment



**Nuclear weapons testing**



**Nuclear power plants**


Pu isotopes,  
 $^{237}\text{Np}$ ,  $^{99}\text{Tc}$ ,  
 $^{236}\text{U}$



**Nuclear reprocessing plants**



**Nuclear accidents**



**Nuclear medicine**

# Sources of Pu, Np, U and Tc in the environment

Sources of Pu and Np in the environment						
Source term	$^{238}\text{Pu}$ , Bq	$^{239}\text{Pu}$ , Bq	$^{240}\text{Pu}$ , Bq	$^{241}\text{Pu}$ , Bq	$^{237}\text{Np}$ , Bq	$^{240}\text{Pu}/^{239}\text{Pu}$ atomic ratio
Nuclear weapons testing	$3.3 \times 10^{14}$	$7.4 \times 10^{15}$	$5.2 \times 10^{15}$	$1.7 \times 10^{17}$	$3.9 \times 10^{13}$	~0.19
Burn up of SNAP-9A	$6.3 \times 10^{14}$	-	-	-	-	-
Thule, Greenland, 1968	-	$1 \times 10^{13}$		-	-	-
Palomares, Spain	-	$5.5 \times 10^{10}$		-	-	-
Chernobyl, 1986	$3.0 \times 10^{13}$	$2.6 \times 10^{13}$	$3.7 \times 10^{13}$	$5.5 \times 10^{15}$	-	~0.39
Sellafield reprocessing plant	$1.2 \times 10^{14}$	$6.1 \times 10^{14}$		$2.2 \times 10^{16}$	-	-
La Hague reprocessing plant	$2.7 \times 10^{12}$	$3.4 \times 10^{12}$		$1.2 \times 10^{14}$	-	-

# Sources of Pu, Np, U and Tc in the environment

Sources of $^{99}\text{Tc}$ and $^{236}\text{U}$ in the environment	
Source term	$^{99}\text{Tc}$ released, Bq
Sellafield nuclear reprocessing plant	$1.72 \times 10^{15}$
La Hague nuclear reprocessing plant	$1.54 \times 10^{14}$
Global weapons fallout (1940s-1970s)	$1.40 \times 10^{14}$
Nuclear accident in Chernobyl	$7.5 \times 10^{11}$
Estimated nuclear accident in Fukushima	$>2.5 \times 10^{11}$
Estimated medical application ( $^{99}\text{Mo}$ - $^{99\text{m}}\text{Tc}$ generator)	$<2 \times 10^{10}$
Estimated nuclear power plants	$<1 \times 10^{10}$
Source term	$^{236}\text{U}$ released, Bq
Natural	$8.4 \times 10^{10}$
Anthropogenic	$2.4 \times 10^{15}$


## Significances of Pu, Np, U and Tc determination



- 1) Environmental risk assessment and monitoring
- 2) Nuclear emergency preparedness
- 3) Routine occupational health monitoring
- 4) Nuclear Forensics
- 5) Nuclear decommissioning and waste disposal
- 6) Radioecology and tracer studies


# Distribution characters of Pu, Np, U and Tc in environmental and biological samples

- 1) Levels are very **low** and **vary** with location or sample type
- 2) Often coexist with **matrix elements** (Ca, Mg, Al, V, Ru, Mo...) and **other interfering radionuclides** (Th, Am, Cm...)





# Traditional analytical methods for Pu and Np



## Our objectives

1. Rapid determination of Pu, Np, U and Tc
2. Automation of the analytical procedure

Specific focuses:

i. Chemical purification


- Protocol simplification and optimization
- Automation

ii. Detection

- Mass spectrometry (ICP-MS, AMS)


# Automation techniques

- Vacuum box
- HPLC
- Flow injection/Sequential injection



Qiao, J. X., Hou, X. L., Miró, M., Roos, P. *Analytica Chimica Acta*. 2009, 652, 66-84.

# Final strategies



# Methods Development-Pu and Np



## Environmental Samples:

- 0.5-200 g of soil, sediment or seaweed
- 50-200 L of seawater

## Parameters optimized:


- Resin type (TEVA, AG1, AG MP-1M)
- Column size (1-20 mL)
- Washing solution (1-8 M HNO<sub>3</sub>)
- Elution solution (NH<sub>2</sub>OHHCl-HCl, 0.1-1.0 M HCl)
- Flow rate (1-5 mL/min)

## **Performance evaluation :**


Chemical yields; <sup>237</sup>Np/<sup>242</sup>Pu chemical yield;

Decontamination of U; Method reability; Sample throughput

# Auto-uint no.1---Sequential injection



Automatically handle 9 samples!  
Work overnight !



Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Analytica Chimica Acta. 2011, 685, 111-119.

# Selected results for soil analysis



Method	Analyte	Resin	Chemical yield of $^{242}\text{Pu}$ , $Y_{\text{Pu}}$ (%)	Chemical yield of $^{237}\text{Np}$ , $Y_{\text{Np}}$ (%)	Ration of $Y_{\text{Np}}/Y_{\text{Pu}}$	$^{239}\text{Pu}$ measured (Bq/kg) *	$^{240}\text{Pu}$ measured (Bg/kg)**	Decontamination factor ***		
								$^{238}\text{U}$	$^{232}\text{Th}$	$^{208}\text{Pb}$
Extraction chromatography <sup>1, 2)</sup>	Pu	TEVA (2mL, $0.7 \times 5$ cm)	97.7 ± 3.4	-	-	0.14 ± 0.01	0.09 ± 0.01	$7.5 \times 10^4$	$2.5 \times 10^4$	$1.3 \times 10^5$
	Np & Pu	TEVA (2mL, $0.7 \times 5$ cm)	88.1 ± 3.4	85.7 ± 3.9	0.97	0.14 ± 0.01	0.09 ± 0.01	$1.0 \times 10^4$	$7.0 \times 10^3$	$1.0 \times 10^4$
Anion chromatography <sup>3, 4)</sup>	Pu (&Np)	AG 1-X4 (50-100mesh), (2mL, $0.5 \times 10$ cm)	103.0 ± 5.2	84.8 ± 5.3	0.75	0.14 ± 0.02	0.09 ± 0.01	$3.9 \times 10^3$	$2.4 \times 10^4$	$2.7 \times 10^4$
		AG 1-X4 (100-200mesh), (2mL, $0.5 \times 10$ cm)	91.6 ± 4.6	75.8 ± 4.6	0.77	0.14 ± 0.01	0.10 ± 0.01	$6.9 \times 10^3$	$1.7 \times 10^4$	$1.0 \times 10^3$
	Np & Pu	AG MP-1M (100-200mesh), (2mL, $0.5 \times 10$ cm)	86.5 ± 4.3	85.3 ± 4.3	0.99	0.14 ± 0.02	0.10 ± 0.01	$3.9 \times 10^3$	$2.5 \times 10^5$	$1.0 \times 10^3$

10 g of soil was used in each analysis. \*The reference value is  $0.140 \pm 0.008$  Bg/kg. \*\*The reference value is  $0.098 \pm 0.006$  Bg/kg. \*\*\* The relative standard deviations were in all instances better than 10%.


1) Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Analytical Chemistry. 2009, 81, 8185-8192.

2) Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Journal of Analytical Atomic Spectrometry. 2010, 25, 1769-1779.

3) Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Analytica Chimica Acta. 2009, 652, 66-84.

4) Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Talanta. 2011, 84, 494-500.

# Optimized chemical purification for Pu and Pu simultaneous determination



Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Analytical Chemistry. 2011, 83, 374-381.



## Performance comparison

Item	Extraction chromatography	Anion exchange chromatography
Price of resin	☹ High (e.g. 5600 €/500 g)	😊 Relatively Low (e.g. 330-1000 €/500 g)
Chemical yields	80-100% (Pu), ☹ 40-80% (Np&Pu)	😊 80-100% (Pu), 70-90% (Np&Pu, AG MP-1M resin)
Separation time	😊 1.5 hr/sample	☹ 2.5-3.5 hr/sample
Decontamination	😊 High ( $1-10 \times 10^4$ for $^{238}\text{U}$ )	☹ Medium ( $1-10 \times 10^3$ for $^{238}\text{U}$ )
Accuracy	😊 High ( $\text{RSD} \leq 5\%$ )	☹ Medium ( $\text{RSD} \leq 10\%$ )
Consumption of chemicals	😊 Low (e.g. 10 mL of conc. $\text{HNO}_3$ /sample)	☹ High (e.g. 80 mL of conc. $\text{HNO}_3$ /sample)
<b>Recommendation</b>	<b>Pu determination using TEVA resin</b>	<b>Pu &amp;Np simultaneous determination using AG MP-1M resin</b>

Qiao, J. X. Risø-PhD-75 (EN), ISBN 978-87-550-3889-9, March 2011, Roskilde, Denmark.

# Methods Development-Pu and Np



## Biological samples:

- 1-5 L urine

## Parameters optimized:


- Co-precipitation techniques ( $\text{Ca}_3(\text{PO}_4)_3$ ,  $\text{Fe}(\text{OH})_3$ ,  $\text{MnO}_2$ , etc.)
- Decomposition of organic matter (acid digestion/dry ashing)
- Washing solution (0.2-1 M  $\text{HNO}_3$ )
- Elution solution (0.025-0.5 M  $\text{HCl}$ )

## **Performance evaluation :**

Chemical yields;  $^{237}\text{Np}/^{242}\text{Pu}$  chemical yield ratio;


Method reability; Sample throughput

# Auto-unit no.2---LOV bead injection



Qiao, J. X., Hou, X. L., Roos, P., Miró, M. Analytical Chemistry. 2013, 85, 2853-2859.

# Analytical procedure for urine analysis



# Selected results for urine analysis



Group no.	Pre-concentration method	Organic matter decomposition	Valence adjustment reagents	Operation time	Chemical yield		
					$^{242}\text{Pu}, \%$	$^{237}\text{Np}, \%$	$^{237}\text{Np}/^{242}\text{Pu}$
1	$\text{Ca}_3(\text{PO}_4)_2$ co-precipitation	Dry ash	Ascorbic acid / conc. $\text{HNO}_3$	13 hr	$84.7 \pm 5.7$	$80.9 \pm 10.7$	0.95
		Acid digestion		8 hr	$46.8 \pm 4.1$	$8.3 \pm 5.4$	0.18
2	$\text{Fe}(\text{OH})_2/\text{Fe}(\text{OH})_3$ co-precipitation	Dry ash	$\text{Fe}/\text{K}_2\text{S}_2\text{O}_5 / \text{conc.}\text{HNO}_3$	6 d	$84.3 \pm 15.6$	$73.3 \pm 33.0$	0.87
		Acid digestion		5.5 d	$80.3 \pm 9.9$	$77.9 \pm 10.9$	0.97
		Acid digestion		6 hr	$51.3 \pm 0.2$	$57.5 \pm 8.8$	1.12
3	<b><math>\text{MnO}_2</math> co-precipitation</b>	<b>Acid digestion</b>	<b><math>\text{Fe}/\text{K}_2\text{S}_2\text{O}_5 / \text{conc.}\text{HNO}_3</math></b>	<b>6 hr</b>	<b><math>88.4 \pm 8.0</math></b>	<b><math>91.4 \pm 10.0</math></b>	<b>1.03</b>
4	$\text{Ca}(\text{OH})_2/\text{Fe}(\text{OH})_2 / \text{Fe}(\text{OH})_3$ co-precipitation	Acid digestion	Ascorbic acid / conc. $\text{HNO}_3$	6 hr	$87.3 \pm 6.6$	$51.2 \pm 1.6$	0.59
5	Evaporation	Dry ash + acid leaching	$\text{Fe}/\text{K}_2\text{S}_2\text{O}_5 / \text{conc.}\text{HNO}_3$	1.5 d	$75.5 \pm 2.6$	$81.1 \pm 3.6$	1.07

# Methods Development-Pu, Np and U



## Environmental Sample:

- 10 L of seawater

## Parameters optimized:


- Resin type (TEVA, UTEVA)
- Washing solution (1-4 M HNO<sub>3</sub>)
- Decomposition of organic matter

## **Performance evaluation :**

Chemical yields; <sup>237</sup>Np/<sup>242</sup>Pu chemical yield ratio;


Method reability; Sample throughput

## Auto-unit no.3---Dual-column sequential injection



Flexible control the connection of two columns!

Qiao, J. X., Hou, X. L., Steier, P., Golser, R. Analytical Chemistry. 2013, submitted.





## Selected results for Pu/Np/U seawater analysis

### Typical analytical performance

Analytical time, h	Chemical yield, %				DU
	$^{242}\text{Pu}$	$^{237}\text{Np}$	$^{237}\text{Np}/^{242}\text{Pu}$ ratio	$^{238}\text{U}$	
8	$73.6 \pm 9.8$	$73.9 \pm 5.6$	$1.01 \pm 0.21$	$97.6 \pm 20.6$	$(5.3 \pm 0.5) \times 10^4$

### Method application with the use of AMS measurement

Sample ID	$^{236}\text{U}/^{238}\text{U}$ , $\times 10^{-8}$	$^{238}\text{U}$ , $\mu\text{g/L}$	$^{236}\text{U}$ atom/L,	Measured value, mBq/L		Expected value, mBq/L	
				$^{237}\text{Np}$	$^{239}\text{Pu}$	$^{237}\text{Np}$	$^{239}\text{Pu}$
North Atlantic-1	$8.88 \pm 1.33$	$2.76 \pm 0.41$	$(6.21 \pm 0.93) \times 10^8$	<0.001	<0.005	-	-
North Atlantic-2	$2.03 \pm 0.30$	$2.17 \pm 0.33$	$(1.11 \pm 0.17) \times 10^8$	0.18	1.02	0.20	1.12
Roskilde Fjord-1	$1.40 \pm 0.21$	$1.65 \pm 0.28$	$(6.88 \pm 1.03) \times 10^7$	<0.001	0.03	-	-
Roskilde Fjord-2	$1.65 \pm 0.25$	$1.65 \pm 0.28$	$(5.85 \pm 0.88) \times 10^7$	0.16	1.18	0.15	1.12

# Methods Development-Pu, Np, U and Tc



## Environmental Sample:

- 200 L of seawater

## Parameters optimized:


- Resin type (TEVA, UTEVA, AG MP-1M)
- Selection of Redox reagents
- Decontamination of interferences

## **Performance evaluation :**


Chemical yields;  $^{237}\text{Np}/^{242}\text{Pu}$  chemical yield ratio;

Method reability; Sample throughput

# Auto-uint no.4---Flow injection



Simultaneously handle 4 samples!



# Summary

Objectives	Achievement
1. Rapid determination of Pu, Np, U and Tc	<ul style="list-style-type: none"> <li>•Environmental solids: 2-5 h/sample</li> <li>•Large volume seawater: 1-2 days/sample</li> <li>•Biological samples: 6 h/sample</li> </ul>
2. Automation of the analytical procedure	<ul style="list-style-type: none"> <li>•Sample pre-concentration: batchwise</li> <li>•Chemical purification: automated                             <ol style="list-style-type: none"> <li>1. Auto-unit no.1: sequentially 9 samples</li> <li>2. Auto-unit no.2: automated column packing</li> <li>3. Auto-unit no.3: automated dual connection</li> <li>4. Auto-unit no.4: simultaneously handle 4 samples</li> </ol> </li> <li>•Measurement: automated</li> </ul>

# Conclusions and perspectives



## Innovation of the previous work:

- Automatic
- Rapid and simple
- No need of Np isotopic tracer
- Low consumption of resins
- High sample throughput
- Low labor intensity

## On-going projects:

- Tracer application studies of Pu and  $^{236}\text{U}$
- Multi-radionuclide determination (Pu/Np,U, Th, Am) in environmental samples

# Thank you!

Jixin Qiao

[jiqui@dtu.dk](mailto:jiqui@dtu.dk)