



Tem holder for sample transfer under reaction conditions

Damsgaard, Christian Danvad; Zandbergen, Henny W.; Wagner, Jakob Birkedal

Publication date:
2013

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):
Damsgaard, C. D., Zandbergen, H. W., & Wagner, J. B. (2013). Tem holder for sample transfer under reaction conditions. Abstract from PICO 2013, Juelich, Germany.

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

TEM HOLDER FOR SAMPLE TRANSFER UNDER REACTION CONDITIONS

Christian D. Damsgaard^{§†}, Henny W. Zandbergen[‡], and Jakob B. Wagner[§]

[§]Center for Electron Nanoscopy, Technical University of Denmark, DK-2800 Lyngby, Denmark, christian.damsgaard@cen.dtu.dk

[†]CINF, Department of Physics, Technical University of Denmark, DK-2800 Lyngby, Denmark

[‡]Kavli Institute of Nanoscience, Delft University of Technology, Lorentzweg 1, 2628 CJ Delft, The Netherlands

Abstract

Environmental transmission electron microscopy (ETEM) studies are usually performed using conventional sample holders in a dedicated ETEM [1] or in a traditional TEM by use of a dedicated high-pressure cell sample holder [2]. In both cases, the setup defines the conditions regarding gas, pressure and temperature, which usually are far from the operando conditions of e.g. heterogeneous catalysis. Our efforts focus on bridging these gaps by establishing *in situ* sample transfer between complementary measurement techniques.

To fully exploit the capabilities of ETEM complementary experiments and characterization techniques are beneficial. Normally, the complementary measurements are done in parallel with experiments separated in time and space [3] or by mimicking a reactor bed by changing the feed gas composition according to reactivity and conversion measured in dedicated catalyst set-ups [4]. Furthermore, dedicated transfer holders have been used to transfer catalyst samples between reactor set-ups and TEM at room temperature in inert atmosphere [5]. To take the full advantage of complementary *in situ* techniques, transfer under reactions conditions is essential.

This study introduces the *in situ* transfer concept by use of a dedicated TEM transfer holder capable of enclosing the sample in a gaseous environment at temperatures up to approx. 900°C. By oxidation and reduction experiments of Cu nanoparticles it is shown possible to keep the reaction conditions during transfer outside the microscope.

This opens up for the possibility to do complimentary *in situ* experiments of the exact same sample without changing the sample condition during transfer. As an example a commercially used methanol catalyst Cu/Zn/Al₂O₃ [6], was reduced at 1 bar in 2% H₂ in He at 220°C and in an *in situ* x-ray diffractometer (XRD) setup and subsequently successfully transferred under reaction conditions to the ETEM.

1. T. W. Hansen et al, Science 294 (2001) 1508
2. J. Creemer et al., Ultramicroscopy 108 (2008) 993
3. P.L. Hansen et al., Science 295 (2002) 2053
4. S. Chenna et al., ChemCatChem 3 (2011) 1051
5. P.J. Kooyman et al., Catalysis letters 74 (2001) 49
6. C. Bates et al. Journal of Catalysis 258 (2008) 334
7. We kindly acknowledge the support from the Danish Ministry of Science and Innovation through the Catalysis for Sustainable Energy Initiative (CASE) at DTU and by The Danish National Research Foundation through CINF DTU. Furthermore, the A. P. Møller and Chastine Mc-Kinney Møller Foundation is gratefully acknowledged for its contribution towards the establishment of DTU Cen.